





# B. ALI IMPEX

Khadim Ali Road, Kotli Behram,  
Sialkot-51310, Pakistan  
Tel: 052-4290628 Fax:052-4290737  
email: info@baliimpex.com  
baliimpex.com



A sign of quality





Cat#	Description	Picture
542/SAN	Diagnostic Set Screw fitting including: 1 Otoscope head double lens with 3 reusable speculas 1 Ophthalmoscope head with lens selection wheel (0 to 20 and 0 to -20 diopters single large circle 1 bent arm illuminator 1 spreadable nasal specula with spring action 1 tongue blade holder 2 laryngeal mirrors 1 battery handle 1 spare bulb	
542/SANL	Same as above but with LED lamp having white bright light	See above
542/ER	Diagnostic Set including: 1 Otoscope head with 3 reusable speculas 1 Ophthalmoscope head with lens selection wheel (0 to 40 and 0 to -25 diopters single large circle 1 bent arm illuminator 1 nasal specula 1 tongue blade holder 1 tongue blade 2 laryngeal mirrors 1 battery handle 1 adaptor for single use speculas 8 single use speculas	
542/ERL	Same as above but with LED lamp having white bright light	See above
542/A	Diagnostic Set, Pin type including: 1 Otoscope head with 3 reusable Speculas. 1 Ophthalmoscope head with lens selection wheel (+0 to +20 and -0 to -20) diopters single large circle, aperture, 1 Bent Arm illuminator, 1 Nasal Speculum 1 Tongue blade holder, 1 tongue Blade 2 Laryngeal mirrors & 1 Battery handle	
542/AL	Same as above but with LED lamp have white bright light	See above
543/FR	Diagnostic Set including: 1 Otoscope head with 3 reusable speculas 1 Ophthalmoscope head with lens selection wheel (0 to 40 and 0 to -25 diopters single large circle 1 bent arm illuminator 1 nasal specula 1 tongue blade holder 1 tongue blade 2 laryngeal mirrors 1 battery handle	
543/FL	Same as above but with LED lamp have white bright light	See above

# B. ALI IMPEX

Khadim Ali Road, Kotli Behram,  
Sialkot-51310, Pakistan  
Tel: 052-4290628 Fax:052-4290737  
email: info@baliimpex.com  
baliimpex.com



A sign of quality

Cat#	Description	Picture
543/GR	<p>Diagnostic Set including:</p> <ul style="list-style-type: none"> <li>1 Otoscope head with 3 reusable speculas</li> <li>1 Ophthalmoscope head with lens selection wheel (0 to 40 and 0 to -25 diopters single large circle)</li> <li>1 bent arm illuminator</li> <li>1 nasal specula</li> <li>1 tongue blade holder</li> <li>1 tongue blade</li> <li>2 laryngeal mirrors</li> <li>1 battery handle</li> </ul>	
543/GL	Same as above but with LED lamp have white bright light	See above
542/G	<p>Otoscope/ ophthalmoscope Set including:</p> <ul style="list-style-type: none"> <li>1 Otoscope head with 3 reusable speculas</li> <li>1 Ophthalmoscope head with lens selection wheel (0 to 40 and 0 to -25 diopters single large circle)</li> <li>1 battery handle</li> </ul>	
542/GL	Same as above but with LED lamp have white bright light	See above
542/J	<p><b>Mini Oto-Ophthalmoscope Set:</b></p> <p>1 Fiber Optic Otoscope head with 10 disposable Speculas. 1 Ophthalmoscope head with lens selection wheel (+0 to +20 and -0 to -20) diopters single large circle.</p> <p>Packed in plastic box</p>	
542/JL	Same as above but with LED lamp have white bright light	See above
543/HR	<p>Otoscope/ ophthalmoscope Set including:</p> <ul style="list-style-type: none"> <li>1 Fiber Optic Otoscope head with 8 reusable speculas</li> <li>1 Spare Lamp</li> <li>1 battery handle</li> </ul>	
542/HL	Same as above but with LED lamp have white bright light	See above

# B. ALI IMPEX

Khadim Ali Road, Kotli Behram,  
Sialkot-51310, Pakistan  
Tel: 052-4290628 Fax:052-4290737  
email: info@baliimpex.com  
baliimpex.com



A sign of quality


Cat#	Description	Picture
542/K	<p><b>Veterinary Oto- ophthalmoscope set:</b> Lock fitting including 1 Otoscope head with 3 reuseable Speculas. 1 Ophthalmoscope head with lens selection wheel (+0 to +20 and -0 to -20) diopters single large circle. Packed in plastic box</p>	
542/KL	Same as above but with LED lamp have white bright light	See above
543/JR	<p>Rechargeable Diagnostic Set including: 1 Otoscope head with 3 reusable speculas 1 Ophthalmoscope head with lens selection wheel (0 to 40 and 0 to -25 diopters single large circle 1 bent arm illuminator 1 nasal specula 1 tongue blade holder 1 tong blade 2 laryngeal mirrors 1 Rechargeable battery handle 1 Charger</p>	
543/JL	Same as above but with LED lamp have white bright light	See above
542/L	<p><b>Oto- ophthalmoscope set:</b> Lock fitting including 1 Otoscope head with 3 reuseable Speculas. 1 Ophthalmoscope head with lens selection wheel (+0 to +20 and -0 to -20) diopters single large circle. Packed in plastic box</p>	
542/LL	Same as above but with LED lamp have white bright light	See above
543/KR	<p>Diagnostic Set including: 1 Fiber Otoscope head with 14 disposable speculas 1 bent arm illuminator 1 tongue blade holder 1 tong blade 2 laryngeal mirrors 1 battery handle</p>	
543/KL	Same as above but with LED lamp have white bright light	See above

# B. ALI IMPEX

Khadim Ali Road, Kotli Behram,  
Sialkot-51310, Pakistan  
Tel: 052-4290628 Fax:052-4290737  
email: info@baliimpex.com  
baliimpex.com



A sign of quality

Cat#	Description	Picture
#542/HN	Mini Otoscope/ Mini ophthalmoscope Set including: 1 Otoscope head fiber optic with 10 disposable speculas 1 Ophthalmoscope head with lens selection wheel (0 to 40 and 0 to -25 diopters 5 different apertures, green filter)  1 battery handle	
542/HNL	Same as above but with LED lamp have white bright light	See above