



### Fowler Bed (BAI-786-501)

Three cranks manual hospital bed  
Dimension:(78x36x20)  
Made of steel powder coated frame  
Bed surface ABS Mesh  
Detachable ABS head and foot board  
ABS side rail for protection  
3 sets manual crank mechanism  
Bumpers on four corners  
2 sets drainage bag holder at both side of the bed  
4 pcs luxurious wheels, with brake



### Fowler Bed (BAI-786-502)

Dimension; 2080\*960\*500mm  
Bedding area covered with M.S. Strips Sheet  
Head & Foot Board ABS, Detachable,  
ABS side railing  
Caster & Locking 125mm, individual brake



### Fowler Bed (BAI-786-503)

Dimension ; Dimensions: (78 x 36 x 24) inches  
Mobile on 5" diagonally locking castors  
Available in strip and wire mesh platform  
Removable head & foot panels  
Powder coated



## Fowler Bed (BAI-786-504)

Full Fowler Bed  
Two Crank Operation  
Removable head and foot panels  
Rotary buffers for protection  
Mobile on 5" diagonally locking castors  
Dimensions: (78 x 36 x 22) inches  
Epoxy powder coated



## Fowler Bed (BAI-786-505)

Semi Fowler Bed  
Single Crank Operation  
Removable Head and Foot panels  
Mobile on 5" diagonally locking castors



## Bedside Locker (BAI-786-506)

Dimension; 18" W x 15" D x 34" H  
Powder coated sheet steel construction  
drawers with one bottom compartment  
Plastic top laminated with raised edges on 3 sides.  
75mm twin castors rust proof waterproof, two lockable



## Bedside Cabinet (BAI-786-507)

Dimension 18" W x 15" D x 34" H  
Top stainless steel with raised edges on three sides  
Open shelf between Drawer and locker  
Epoxy powder coated



## Bedside Locker (BAI-786-508)

Dimension 16" L x 16" W x 36" H  
U shape Comprises a drawer and compartment fitted in U-shaped legs,  
Table top fitted with SS sheet.  
Powder coated finish



## Over Bed Table (BAI-786-509)

Top size: 80 x 32 cm

Table top made of Laminated sheet

Fitted on steel tubular telescopic

height adjustment of two rectangular tubular having height adjustment with tightening of screw knob.

Mobile on 50 mm Dia wheels



## Instruments Mayo Table (BAI-786-511)

Stainless steel construction foot pedal-operated lock adjusts height from 36" to 62"

Tray is removable for easy cleaning

2" double ball bearing casters



## Delivery Table (BAI-786-512)

Fully Stainless-Steel Construction argon welded non-foldable sturdy, stable structure

Collapsible side rails

Table top cushioned with high density mattress cover with rexene. Thickness 3 inch (Approx)

Dimensions: 1700 X 500 X 800 mm (Approx)



### Ward Screen (BAI-786-513)

Frame in four (03) parts constructed from fully plastic coated, MS hollow tubular pipe

All frames fixed together with ABS universal joint for easy movement of each part.

Polyester powder coating

Coaster on 75mm



### Mayo Stand (BAI-786-514)

Frame in four (04) parts

MS hollow tubular pipe

Curtain fitted with spring loaded rods

Mobile on six, coaster on 75mm dia

Polyester powder coating



### Instrument Trolley (BAI-786-15)

Stainless Steel.

Dimension 24" x 36"

Mobile on castors



## Instrument Cabinet (BAI-516)

Material SS

Size: Standard Size (~1800mm x900mm x 400 mm)

Double Door see through glass



## Dressing Trolley (BAI-786-517)

Stainless Steel trolley with bucket and bowl attached

Dimension 75 l x 45w x 90 h cm

Wheels: Mobile on castors.



## Instrument Trolley With Drawer (BAI-786-518)

Stainless Steel with two draws

Wheels: 4 wheels. Mobile on castors



## Baby Cot (BAI-786-519)

With mattress size 2 inch

12 degrees Trendelenburg, 6 degrees reverse Trendelenburg. Four 75mm individually braking castors

Powder coated



## **Baby Bed (BAI-786-520)**

Size: 900X600X750mm

Structure made of Mild Steel tube pipe 1"X1" 16 SWG bended not welding

One side fixed and other side should be movable. 3" mattress and 4" wheel with lockable



## **Crash Cart (BAI-786-521)**

DIMENSIONS:26"x26"x59"

Full extensions drawers offer access to entire contents of drawers

Heavy duty construction

Extra top shelf for monitors

Stainless steel adjustable IV Pole

Dual Oxygen cylinders

The Trolley is protected by wrap-around-bumper drawers are furnished with a central locking device for protection

CPR Board included

Epoxy powder coated finish

Mobile on 5" swivel castors





## Patient Stool (BAI-786-522)

Stainless steel structure Rust free 12" dia 18 SWG

Height Adjustable for 43 to 55 cm



## Examination Couch (BAI-786-525)

- Dimension; 180 L x 60 W x 46 H CMS
- Patients Examination
- Couch made of stainless-steel structure
- Manually operated
- Upholstery with Washable material
- Manual Head raise
-



## Examination Couch (BAI-786-526)

Dimension; 180 L x 60 W x 46 H CMS

Patients Examination couch made of Mild steel

Manually operated. Upholstery with washable material

Manual Head raise

Powder coated



## Stretcher Trolley (BAI-786-527)

Dimensions: 75" x 24" x 31" (L x W x H)

Mild Steel construction

Utility basket underneath

Removable stretcher top

*Cylinder holder*

IV Pole supplied

Mobile on 6" castors



## Drip Stand (BAI-786-528)

Heavy Base Five Wheels Movable

4 Sides Hangers SS

Height Adjustable



## Clean Linen Trolley (BAI-786-529)

24 x 48 x 36 Inch

Stainless steel pipe frame

Dual wash able parachute bag easily removable

Mobile on four caster



## Linen Trolley (BAI-786-530)

- Mounted on 100mm non rusting castors two with brakes
- Three shelves
- Linen trolley with horizontal supports in tube with bag
- Mounted on 100mm non rusting castors two with brakes
- Linen trolley with horizontal supports in tube with bag



## Scrub Double Bay (BAI-786-531)

Stainless steel construction of non-magnetic material 304

Double bay system, bay size 5-6'x3'x2'.

Mechanical water valve Stopped Basin anti splash ends

Elbow operating water flow

Soap dispensing wall mounted



## Scrub Unite (BAI-786-532)

Single bay system, bay size 3'x3'x2'.

Stainless steel construction of non-magnetic material 304

Elbow operating water flow



## Wheel Chair (BAI-786-533)

Overall size: 32 x 28 x 20 Inch

Mild steel Construction

Foot rest foldable

Wire mesh utility basket underneath