








Fiber Optic Laryngoscope Sets

Cat#	Description	Picture
L3-150F-S1	<p>Reusable Laryngoscope Set Adult & Ped, MacIntosh and Miller Fiber Optic. Consist of 5 stainless steel Mac Blades(0, 1,2,3 & 4) 3 stainless steel Mil blade (00, 0 & 1) 1 C-size Battery Handle chrome plated 1 AA-size batter handle chrome plated Green Line conform to ISO 7376</p> 	 <p>Brand: DoctorNice</p>
Cat#	Description	Picture
L3-150F-S2	<p>Reusable Laryngoscope Set, MacIntosh Fiber Optic. Consist of 4 stainless steel Mac Blades(1,2,3,4) 1 C-size Battery Handle chrome plated (HF1) Stored in plastic box : Green Line conform to ISO 7376</p>  <p>Brand: DoctorNice</p>	
L3-150F-S3	<p>Reusable Laryngoscope Set, MacIntosh Fiber Optic. Consist of 4 stainless steel Mac Blades(1,2,3,4) 1 C-size Battery Handle chrome plated (HF1) Stored in plastic box : Green Line conform to ISO 7376</p> <p>Brand: DoctorNice</p>	
Cat#	Description	Picture
L3-150F-S4	<p>Reusable Laryngoscope Set, MacIntosh Fiber Optic. Consist of 3 stainless steel Mac Blades(2,3,4) 1 C-size Battery Handle chrome plated(HF1) Stored in plastic box : Green Line conform to ISO 7376</p> 	 <p>Brand: DoctorNice</p>

Fiber Optic Laryngoscope Sets

Cat#	Description	Picture
L3-155F-S1	<p>Reusable Laryngoscope Set, Miller Fiber Optic. Consist of 5 stainless steel Miller Blades(0,1,2, 3 & 4) 1 C-size Battery Handle chrome plated(HF2) Stored in plastic box: Green Line conform to ISO 7376</p> 	 <p>Brand: </p>
L3-155F-S2	<p>Reusable Laryngoscope Set, Miller Fiber Optic. Consist of 3 stainless steel Miller Blades(0,1,2) 1 AA-size Battery Handle chrome plated(HF2) Stored in plastic box: Green Line conform to ISO 7376</p>  <p>Brand: </p>	<p>Picture</p> 
L3-153F-S2	<p>Reusable Mega Fiber Optic Laryngoscope Set consists of : 4 Mega Laryngoscope Blade (L3-153F) size 2, 3, 4, 5)La 1 Laryngoscope Battery Handle Stored in Plastic Box: Green Line conform to ISO 7376</p>  <p>Brand: </p>	
L3-154F-S1	<p>Reusable Mc Coy Fiber Optic Laryngoscope Set consists of : 3 Mc Coy Fiber optic Laryngoscope Blades (L3-154F) sizes 2, 3 & 4 1 Laryngoscope Battery Handle, C-Size (HF1) Green Line conform to ISO 7376</p>  <p>Brand: </p>	<p>Picture</p> 


Fiber Optic Laryngoscope Sets

Cat#	Description	Picture
L3-154F-S2	<p>Reusable Mc Coy Fiber Optic Laryngoscope Set consists of : 2 Mc Coy Fiber optic Laryngoscope Blades (L3-154F) sizes 3 & 4 1 Laryngoscope Battery Handle, C-Size (HF1) Green Line conform to ISO 7376</p>  <p>Brand: </p>	
L3-154F-S3	<p>Reusable Mc Coy Fiber Optic Laryngoscope Set consists of : 1 Mc Coy Fiber optic Laryngoscope Blades (L3-154F) sizes 4 1 Laryngoscope Battery Handle, C-Size (HF1) Green Line conform to ISO 7376</p>  <p>Brand: </p>	
L3-164F-S1	<p>Mac Gold line Laryngoscope Set, Fiber Optic consists of: 3 Blades (L3-164F) size 2, 3, 4, Fiber optic light carrier can be easily detached from the blade and attached to the blade which facilitate cleaning of the blades 1 Water Proof Battery Handle (HF7), C-size, brass chrome plated Store in plastic box .</p>  <p>Brand: </p>	<p>Picture of set is not available but pictures of blades & handles are as under:</p> 
L3-164F-S2	<p>Miller Gold line Laryngoscope Set, Fiber Optic consists of: 3 Blades (L3-164F) size 2, 3, 4, Fiber optic light carrier can be easily detached from the blade and attached to the blade which facilitate cleaning of the blades 1 Water Proof Battery Handle (HF7), C-size, brass chrome plated Store in plastic box .</p>  <p>Brand: </p>	<p>Picture of set is not available but pictures of blades & handles are as under:</p> 

Conventional Laryngoscope Sets

Cat#	Description	Picture
L3-150-S1	<p>Reusable Laryngoscope Set, MacIntosh Conventional. Consist of 5 stainless steel Mac Blades(0,1,2,3,4) 1 C-size Battery Handle chrome plated Stored in plastic box: Conform to ISO 7376</p>  <p>Brand: Doctor Nice</p>	
L3-150-S2	<p>LED Reusable Laryngoscope Set, MacIntosh Conventional. Consist of 4 stainless steel Mac Blades(0,1,2,3,4) 1 C-size Battery Handle chrome plated Stored in plastic box : Conform to ISO 7376</p>  <p>Brand: Doctor Nice</p>	
L3-150-S3	<p>Reusable Laryngoscope Set, MacIntosh Conventional. Consist of 3 stainless steel Mac Blades(2,3,4) 1 C-size Battery Handle chrome plated Stored in plastic box : Conform to ISO 7376</p>  <p>Brand: Doctor Nice</p>	

Conventional Laryngoscope Sets

Cat#	Description	Picture
L3-150-S4	<p>Reusable Laryngoscope Set, MacIntosh Conventional. Consist of 3 stainless steel Mac Blades(2,3,4) 1 spare lamp 1 C-size Battery Handle chrome plated Stored in artificial leather pouch Conform to ISO 7376 Brand: </p>	
L3-155-S1	<p>Reusable Laryngoscope Set, Miller Conventional. Consist of 4 stainless steel Mac Blades(0,1,2,3) 1 AA-size Battery Handle chrome plated Stored in plastic box: Conform to ISO 7376  Brand: </p>	
L3-155-S2	<p>Reusable Laryngoscope Set, Miller Conventional. Consist of 3 stainless steel Mac Blades(00,0,1) 1 AA-size Battery Handle chrome plated Stored in plastic box : Conform to ISO 7376  Brand: </p>	
L3-167-S1	<p>Reusable Conventional Mc Coy Laryngoscope Set consists of : 2 Conventional Mc Coy Laryngoscope blades (L3-167) size 3 & 4. 1 Laryngoscope Battery Handle Stored in plastic box: Conform to ISO 7376  Brand: </p>	