

B. ALI IMPEX

Khadim Ali Road, Kotli Behram,
Sialkot-51310, Pakistan
Tel: 052-4290628 Fax:052-4290737
email: info@baliimpex.com web: www.baliimpex.com



A sign of quality

FIBER OPTIC LARYNGOSCOPES

#L3-153F

Reusable Mega Fiber Optic Laryngoscope Blades. Cleanliness of these blades is very easy because of their design (inside fiber optic guide). These blades are made of Non-Magnetic stainless steel, Sizes :2, 3, 4, 5
Green Line conform to ISO 7376



#L3-154F

Reusable Mc Coy Fiber Optic Laryngoscope Blades. The blades are for difficult intubation. These are made of Non-Magnetic stainless steel. Sizes: 3, 4 & 5.
Green Line conform to ISO 7376



#L3-150F

Reusable Professional MacIntosh Fiber Optic Laryngoscope Blade, Made of Grade 304 Stainless Steel (Non-Magnetic), fitted with Glass Fiber Optic light guide covered with stainless steel pipe. Sizes: 0,1, 2, 3,4, 5.
Green Spec



#L3-152FD

Disposable MacIntosh Fiber Optic Laryngoscope Blades. These blades are made of Non-Magnetic stainless steel with aluminum base and fitted with plastic Fiber Optic light guide. Sizes: 0, 1, 2, 3, 4, 5. Green Line conform to ISO 7376




#L3-151F

Single use MacIntosh Fiber Optic Laryngoscope Blades. These blades are made of Grade 304 Stainless Steel (Non-Magnetic) with plastic heal fitted with plastic Fiber Optic light guide. Sizes available 0, 1, 2, 3,4, 5. Green Spec



#L3-167F

Reusable Professional Miller-fix Fiber Optic Laryngoscope Blades. Fiber optic light carrier has been fixed. These blades are made of Grade 304 Stainless Steel (Non-Magnetic), fitted with Glass Fiber Optic light guide covered with stainless steel. Sizes available 1, 2, 3, 4. Green Spec.




#L3-155F

Reusable Professional Miller Fiber Optic Laryngoscope Blades. Fiber optic light carrier has been fixed. These blades are made of Grade 304 Stainless Steel (Non-Magnetic), fitted with Glass Fiber Optic light guide covered with stainless steel. Sizes available 00, 0, 1, 2, 3, 4. **Green Spec.**



#L3-167F

Disposable MacIntosh Fiber Optic Laryngoscope Blades. These blades are made of Non-Magnetic stainless steel with aluminum base and fitted with plastic Fiber Optic light guide. Sizes: 0, 1, 2, 3, 4, 5.
Green Line conform to ISO 7376



L3-156F

Single use Miller Fiber Optic Laryngoscope Blades. These blades are made of Grade 304 Stainless Steel (Non-Magnetic) with plastic heal fitted with plastic Fiber Optic light guide. Sizes available 00, 0, 1, 2, 3, 4
Green Spec




HF1

Battery Handle medium/C-size, made of brass chrome plated fitted with Xenon lamp or LED, Green Spec.



HF2

Battery Handle AA-size, made of brass chrome plated fitted with Xenon lamp or LED. Green Spec




HF5

Battery Handle AA-size, made of brass chrome plated fitted with Xenon lamp or LED. Green Spec



HF0: Disposable Battery Handle C-size, made of Aluminum, Xenon Lamp or LED Green Spec.

HF00: Disposable Battery Handle C-size, made of plastic + aluminum, Xenon Lamp or LED Green Spec.



HF0



HF00

B. ALI IMPEX

Khadim Ali Road, Kotli Behram,
Sialkot-51310, Pakistan
Tel: 052-4290628 Fax:052-4290737
email: info@baliimpex.com web: www.baliimpex.com



A sign of quality

CONVENTIONAL LARYNGOSCOPES

L3-167
Reusable Conventional I Mc Coy Conventional Laryngoscope Blades. The blades are for difficult intubation. These blades are made of Non-Magnetic Stainless steel. **Sizes: 3, 4 & 5 Conform to ISO 7376**



L3-167-LD
Disposable Conventional Mc Coy Conventional Laryngoscope Blades. 2.5V lamp. The blades are for difficult intubation. These blades are made of Non-Magnetic Stainless steel. **Sizes: 3, 4 & 5 Conform to ISO 7376.**



L3-150
Re-usable Conventional Laryngoscope Blade, MacIntosh, made of 304 grade stainless steel (non-magnetic) with 2.5 V lamp, sizes available: 0, 1, 2, 3, 4, 5.




L3-155
Re-usable Conventional Laryngoscope Blade, Miller, made of 304 grade stainless steel (non-magnetic) with 2.5 V lamp, sizes available: 00, 0, 1, 2, 3, 4.



L3-151D
Disposable Conventional Laryngoscope Blade, MacIntosh, made of 304 grade stainless steel with aluminium base, 2.5 V lamp, sizes available: 0, 1, 2, 3, 4, 5.



L3-156D
Disposable Conventional Laryngoscope Blade, Miller, made of 304 grade stainless steel with aluminium base, 2.5 V lamp, sizes available: 00, 0, 1, 2, 3, 4.




L3-170-L
LED Single use Ready Scope Conventional Laryngoscope system (Blade MacIntosh, all stainless steel of grade 304 (non-magnetic) with LED lamp. **These blades are supplied with handle made of aluminium having on/off button. This product is ever ready for use, just push the button to illuminate the blade. Blade sizes: 0,1,2,3,4 & 5. Supplied with battery cells**



L3-170-L2G
2G Mediscope, Conventional, LED Laryngoscope system (Blade MacIntosh, all stainless steel of grade 304 (non-magnetic) with LED lamp. **These blades are supplied with handle made of aluminium. This product is ever ready for use, just pull the blade to "on" position to illuminate the blade. Blade. Sizes available: 0,1,2,3,4 & 5. Supplied with battery cells**




L3-152MH
Conventional Laryngoscope MacIntosh, all stainless steel of grade 304 (non-magnetic) with 2.5 V lamp. **These blades are supplied with AA size handle having on/off button. This product is ever ready for use, just push the button to illuminate the blade. Blade sizes: 0,1,2,3,4 & 5**



L3-157MH
Conventional Laryngoscope Miller, all stainless steel of grade 304 (non-magnetic) with 2.5 V lamp. **These blades are supplied with AA size handle having on/off button. This product is ever ready for use, just push the button to illuminate the blade. Blade sizes:00, 0,1,2,3 & 4**



HC4
Battery Handle, small/AA- size, made of brass, chrome plated for conventional laryngoscope blades
HC5
Battery Handle, Medium/C- size, made of brass, chrome plated for conventional laryngoscope blades



HC8
Battery Handle, Medium/C size, made of brass, chrome plated for conventional laryngoscope blades
HC0
Disposable Battery Handle, Medium/C size, made of Aluminium for conventional laryngoscope blades
HC8
Battery Handle, Stubby, made of brass, chrome plated for conventional laryngoscope blades



SETS OF THE LARYNGOSCOPE CAN BE COMPILED ACCORDING TO YOUR CHOICE (QUANTITY OF BLADES AND HANDLES) IN THE FOLLOWING BOX. BOX PRICE WILL BE US\$ 3.20 EACH.

