



Sign of Quality

# B. Ali Impex

Manufacturers & Exporters



### #560. Fiber Optic mini Oscope

Available in different colours lens x 4 Battery cells size AA(not supplied). On/Off Switch Pneumatic Otoscopy port, 10 single use Speculas.



### #560B: Fiber Optic mini Oscope:

Available in different colours lens x 4 Battery cells size AA(not supplied). On/Off Switch Pneumatic Otoscopy port, 12 single use Specula.



### #560C: Fiber Optic mini Oscope:

Different colours lens x 4, Battery cells size AA(not supplied). On/Off Switch Pneumatic Otoscopy port, 12 single use specula.




### 570/RS

**Conventional Mini Oscope** with metal head. Different colours, Handle without cup, AA Size, with on.off switch, 3 reusable specula, 10 single use speculas and 1 attachment for single use speculas **Soft pouch.**



**# L1-111/A: Oscope** Pin contact/screw fitting, made of brass, chrome plated, lens x 4, metal pneumatic otoscopy port, rheostat button, three black plastic reusable speculas.



**# L1-110/B: Oscope** lock fitting, made of brass, chrome plated, lens x 4, metal pneumatic otoscopy port, rheostat button, three black plastic reusable speculas



**#L1-111: Oscope** Pin contact/ screw fitting, made of brass, chrome plated, double lens (lens x 1 & lens x 4), metal pneumatic otoscopy port, rheostat button, three black plastic reusable speculas.



**#L1-110: Oscope** lock fitting, made of brass, chrome plated, double lens (lens x 1 & lens x 4), metal pneumatic otoscopy port, rheostat button, three black plastic reusable speculas



**#541/C: Oscope** lock fitting, made of brass, chrome plated, lens x 4, metal pneumatic otoscopy port, rheostat button, three black plastic reusable Speculas.



**560R: Fiber Optic mini Oscope Reister Style:** Available in different colours lens x 4 Battery cells size AA(not supplied). On/Off Switch Pneumatic Otoscopy port, 12 single use Specula





Sign of Quality

# B. Ali Impex

Manufacturers & Exporters



**# 540: Fiber Optic Otoscope** lock fitting, made of brass, chrome plated, lens x 4, plastic pneumatic otoscopy port, rheostat button, nine black plastic single speculas in a specula case.



**#540/A: Fiber Optic Otoscope** Pin contact/screw fitting, made of brass, chrome plated, lens x 4, plastic pneumatic otoscopy port, rheostat button, nine black single use speculas in a specula case.



**#L1-113: Ophthalmoscope Set:** with lens selection wheel (+0 to +40 and -0 to -25) diopters single large circle and C-size battery handle with rheostat button.



**#L1-116: Ophthalmoscope Set:** with lens selection wheel (+0 to +20 and -0 to -20) diopters single large circle and C-size battery handle with rheostat button. Packed in artificial leather pouch.



**#L1-118: Ophthalmoscope Set:** with lens selection wheel (+0 to +20 and -0 to -20) diopters single large circle and C-size battery handle with rheostat button. Packed in artificial leather pouch.



**#L1-114: Ophthalmoscope Set:** with lens selection wheel (+0 to +40 and -0 to -25) diopters single large circle and C-size battery handle with rheostat button.



**#L1-120: Ophthalmoscope set:** composed of one threaded handle, AA size diagnostic head, simple to use with all the basic functions such as high quality corrective lens diopters from 0 to +20 and from 0 to -20. Apertures include small, large and semi circle, fixation star and green filter to increase contrast.

